




# $\mu$ Phase<sup>®</sup> & $\mu$ Shape

Modular and compact  
interferometers



# μPhase<sup>®</sup> interferometers

## Measure with highest precision

μPhase<sup>®</sup> interferometers offer objective and precise measurement results of surface and wavefront measurements - quick and reliable.

μPhase<sup>®</sup> interferometers are compact, small and lightweight digital tools which can be used in almost any working environment. These measuring devices are perfectly complemented by the μShape measurement and analysis software to fulfill the highest expectations of quality management.

## Strong arguments for μPhase<sup>®</sup> line of interferometers

- Compact size and modularity enable adaptation to a variety of production and working environments
- Objective digital measurement prevents human errors
- Well-structured and comprehensive software supports both production and lab use

## Measuring without leaving marks behind


The μPhase<sup>®</sup> interferometer systems are used for measuring high precision optical components made of glass, plastic, metal, ceramic and similar. The non-contact measurement method prevents damage to the sample under test and gives the most exact evaluation of the entire surface or wavefront.

## Modular system providing stand alone interferometers and turnkey solutions

TRIOPTICS offers μPhase<sup>®</sup> interferometers as self-contained modular parts as well as pre-defined turnkey solutions.

μPhase<sup>®</sup> customers especially appreciate the space saving and modular concept of the μPhase<sup>®</sup> product line which allows the cost-effective utilization of the instruments. The different parts of the μPhase<sup>®</sup> interferometer line are all compatible and form powerful measurement devices.

μShape interferometer software	μPhase <sup>®</sup> sensors				
	μPhase <sup>®</sup> 3.1		μPhase <sup>®</sup> 3.2		μPhase <sup>®</sup> 3.3
	μPhase <sup>®</sup> turnkey solution (sensor + stand + software)				
	PLANO DOWN	PLANO UP SPHERO UP	VERTICAL/ VERTICAL PRO	UNIVERSAL	ST + R
μLens PLANO and SPHERO objectives					



Different samples examples: molding tool, IOL, ceramic seal surface, Zerodur sphere

# μPhase<sup>®</sup> 3 sensors

## Most flexible interferometer sensors

These highly integrated phase-shifting Twyman-Green interferometer sensors meet the toughest demands for modern quality management. In combination with the measurement and analysis software μShape this high-performance precision measuring instrument provides information about the specimen's surface, wavefront or test objective aberration.



## Advantages of μPhase<sup>®</sup> sensors

- Compact size and modularity enable adaptation to different production and working environments
- Simple and fast adoption to different reflectivities for optimal image contrast adjustment
- Wide field of view alignment mode: Simple and fast alignment of the sample due to a second camera for alignment purposes
- High resolution cameras: 1200x1200 pixels and 600x600 pixels (μPhase<sup>®</sup> 3.1)
- Object plane focusing ability (μPhase<sup>®</sup> 3.3 only)
- High flexibility: Useable in any orientation and different stands
- Measurement accuracy traceable to international standards
- Standard measuring wavelength 632.8 nm; customized versions measuring at wavelengths from 355 nm to 1064 are also available upon request
- Robust, dust-proof housing

	Technical data
Measurement technique	Twyman-Green phase-shifting interferometer, convertible to Fizeau measurement mode
Measurement capability	Measurement of surface topography of reflective surfaces and wavefronts of optical systems in transmission
Laser wavelength	632.8 nm; option: wavelength between 355 nm and 1100 nm
PV repeatability <sup>1)</sup>	$\lambda / 400$ ( $\lambda = 632.8$ nm)
RMS repeatability <sup>1)</sup>	$\lambda / 6500$ ( $\lambda = 632.8$ nm)
Measurement uncertainty	$\lambda / 20$ ( $\lambda = 632.8$ nm), other on request
Camera resolution	μPhase <sup>®</sup> 3.1: 600 x 600 pixel μPhase <sup>®</sup> 3.3: 1200 x 1200 pixel
Digitalization	12 bit
Laser class	with safety mount: class 1; Laser itself: 3B

1) Determined from 100 consecutive measurements over 96% clear aperture with 16 phase averages in a stable environment.


# μLens objectives

The μLens PLANO collimated test objectives and the μLens SPHERO spherical objectives complement the μPhase<sup>®</sup> interferometry systems and increase the flexibility and modularity of the complete system. The μLens PLANO objectives allow for measurements of flat surfaces or optical systems in transmission from 2mm up to 150mm in diameter. The spherical test objectives range μLens SPHERO allows testing of spherical and weak aspherical surfaces, as well as optical systems in transmission. Measurable sample radius (convex and concave) and diameter depend on the specific test lens.

The 100mm PLANO objectives also allow to attach common Fizeau objectives (Fizeau lenses / Transmission spheres).

## Further advantages:

- Existing μPhase<sup>®</sup> systems can be expanded easily and at low cost thanks to the modularity and compatibility of the objective design
- Testing of small samples with radii below 1mm is possible
- High measuring precision through minimum wavefront aberration of the μPhase<sup>®</sup> and μLens objectives
- Field of view correction allows high measurement safety and interferometry with high fringe densities



CX MAX = max. radius of curvature for convex surfaces (R > 0)      CS = Radius of curvatures of the sample  
 CX MIN = min. radius of curvature for convex surfaces (R > 0)  
 CC MAX = max. radius of curvature for concave surfaces (R < 0)  
 CC MIN = min. radius of curvature for concave surfaces (R < 0)

μPhase<sup>®</sup> focusing range for imaging of spherical surfaces


## μLens PLANO overview

Description	Numerical Aperture [mm]	Sample diameter [mm]
μLens PLANO 2	2	0,2 - 2
μPhase <sup>®</sup> base module	5	1-5
μLens PLANO 10	10	0.1-10
μLens PLANO 50	50	10-50
μLens PLANO 100	100	20-100
μLens PLANO 150	150	30-150


Information about objectives for combination with μLens PLANO 100 and 150 available on request.

### μLens SPHERO overview

Description	f-number	NA	$\alpha$ [°]	CXmax [mm]
μLens SPHERO 10 f/0.7	0.7	0.71	90	8.0
μLens SPHERO 10 f/1.0	1	0.50	60	13.0
μLens SPHERO 10 f/1.5	1.5	0.34	40	20.0
μLens SPHERO 10 f/3	3.0	0.17	19	43.0
μLens SPHERO 10 f/5.2	5.2	0.1	11	73.0
μLens SPHERO 10 f/7.6	7.6	0.07	4	108
μLens SPHERO 10 f/12	12	0.04	2	123



Description	f-number	NA	$\alpha$ [°]	CXmax [mm]
μLens SPHERO 50 f/0.7	0.7	0.71	90	26
μLens SPHERO 50 f/1.0	1	0.50	60	45
μLens SPHERO 50 f/1.5	1.5	0.34	40	70
μLens SPHERO 50 f/2.4	2.4	0.21	24	130
μLens SPHERO 50 f/4.1	4.1	0.12	14	225




# μPhase<sup>®</sup> turnkey solutions



Standard  Option

	PLANO DOWN	PLANO UP	SPHERO UP
1 Testing of flat surfaces	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
2 Testing of spherical surfaces			<input checked="" type="checkbox"/>
3 Testing of aspheric surfaces			
4 Testing of wavefronts in transmission	<input type="checkbox"/>		
5 Radius of curvature measurement			relative
6 Low vibration sensitivity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
7 Production use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
8 Quality management use	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
9 R&D department use			
10 Measurement setup	vertical	vertical	vertical
11 Modular / upgradeability	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
12 Stage	table device	table device	table device
Special features			
13 Stand-alone setup (no optical table needed)	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
14 Radius-/position measurement			
15 Automated multiple sample measurement			
16 Usage of CGHs of aspheress			
17 Second moveable platform			



	ST	ST + R	VERTICAL	VERTICAL PRO	UNIVERSAL
1	■	■	■	■	■
2	■	■	■	■	■
3			□	□	□
4	□	□	□	□	□
5	relative	relative, absolute	absolute, automated	absolute, automated	absolute
6					
7	■	■	■	■	□
8	■	■	■	■	
9	■	■	■	■	■
10	vertical	vertical	vertical	vertical	horizontal
11	■	■	■	■	■
12	table device	table device	table device	table device	
13	□	□	□	□	
14	■	■	■	■	■
15			□	■	
16					□
17					□

# μShape Software

The μShape software is used for topography measurement of flat, spherical, cylindrical, toric and aspherical surfaces or wavefronts in production, laboratory and research. Add-on modules enable to adapt the software to custom specific demands. These modules can be added at any time even after the initial purchase.


The μShape software is pre-installed on a state-of-the-art PC or laptop, included with every TRIOPTICS' μPhase<sup>®</sup> interferometer system. Ready-to-use configuration guarantees fast start of work.

With its clear and well structured user interface μShape perfectly deals with the variety of measurement requirements and provides several modules which expand the capabilities of μShape.

μShape works with all Windows<sup>®</sup> versions and is designed for ease-of-use and modular functionality. It controls and displays the measurement results, stores and documents all measurement raw data and ensures maximum transparency and traceability.

## Upgrade laboratory interferometer

The μShape Interferometer Software was originally developed for the μPhase<sup>®</sup> compact interferometers. In combination with the our GenPack the μShape software works with third party interferometers, too.





## Features of the $\mu$ Shape measuring and analysis software (selection)


- Different **user levels** with different access rights
- **Shortcuts** for most used program functions
- Comprehensive context-sensitive **direct help**
- Various **program modes** enable the separate visualization of calibration and measurement processes and its parameter with an integrated live camera image
- Automatic updates of displays and images after every change of analysis parameters or new measurement, allows **re-analyses** of measurements without new measurement
- Easily **pre-configured templates** for a wide range of measuring tasks and analyses; templates contain all parameters and settings, including configuration of screen display
- Graphic windows can be stored in several **graphic formats**
- **Export** of individual parameters or of selected data fields as text, binary or other common file formats (e.g. QED, Zygo XYZ, DigitalSurf) for external processing
- The **measuring results** are presented as parameters or graphically as a cross section, in 2D or 3D plots



- Configurable **measurement reports** show results at a glance
- Several analyses like **Zernike, Seidel, ISO, Legendre, Slope**




- Access control and configuration of **add-on modules** by **dongle**
- Analysis of **aspherical, cylindrical** or **toric** surfaces in spherical or CGH setups
- **External interface** for controlling the interferometer by external programs, e.g. in automated systems
- **MTF** analysis of focal or afocal optical components and systems
- Measuring **multiple apertures** in one shot, e.g. on polishing heads
- Statistical analysis of multiple sub-apertures at the same time (**MultiStat**) including tolerating and pass/fail indication
- **Prism** and **wedge** measurement and analysis
- Consideration of known sample deviations e.g. deviations caused by the optical design (**Sample Normal Data**)
- Analysis of the **tool offset** of lathe machines
- Analysis of **wafer plates**
- **Static fringe** analysis for fast one shot measurements in instable environments



# Interferometry


In interferometry coherent wavefronts are superimposed. The result of this superposition is a fringe pattern, the so-called interferogram. In case of two beam interference each fringe represents a constant phase difference between both waves. Thus the interferogram is a kind of a contour map of the test sample.



Interference principle

The standard design of an interferometer for surface shape testing consists of a collimated coherent light source which is divided by a beam splitter into two beams. The test beam is transformed by a beam shaping optic into a wavefront of nearly the same shape as the sample (commonly flat or spherical). Thus the rays of the test beam intersect the sample under test perpendicularly, are reflected in themselves and embossing the shape errors to the test wavefront. The modified test wavefront is recombined by the beam splitter with the reference beam, reflected at the internal interferometer reference surface, and imaged to the camera sensor. The space of both interferometer arms builds the test cavity. The interferometer measures the optical path difference (OPD) of this cavity for each point independently. Two setups are commonly used for surface and wavefront testing.

## Twyman-Green setup



Twyman-Green principle

A Twyman-Green interferometer is a modified Michelson interferometer. This configuration offers high flexibility, because both interferometer arms can be modified independently of each other. So the intensity of reference and test arm can be easily adjusted to each other in order to get maximum fringe contrast. This is necessary when testing samples with different reflectivity and increases the range of applications enormously. Only a maximum fringe contrast enables a maximum resolution in depth.


As reference surface a surface can be used that is inexpensive and accurately producible independent from the sample size. The adaption to the

sample size is done by conventional beam shaping optics introduced to the test arm. Contrary to the beam shaping optics for Fizeau interferometers (see next) these optics do not require an expensive Fizeau surface as final surface.

Using a flat reference surface also guarantees equal phase shift over entire aperture.

As consequence of this flexibility the interference patterns which can be seen are not caused by the sample errors only but also by the aberrations of the additional optics in both interferometer arms. However, nowadays samples are not anymore evaluated according to its fringe pattern but by a determination of the phase map causing the fringe pattern. During this analysis the aberrations of the additional optics can be easily considered. Finally the software provides an objective digital measurement result.

## Fizeau setup



The most commonly used interferometer for surface testing are Fizeau interferometers. The last surface of the beam shaping optic is the so-called Fizeau surface. It must have the same shape as the sample to be tested (commonly spherical or flat) and is placed concentric into the optical path, so the individual rays intersect perpendicular to the Fizeau surface. The major part of the light passes the Fizeau surface and is reflected at the test surface. The returning light interferes with the part of the light reflected at the Fizeau surface. So the Fizeau surface

acts as beam splitter as well as reference surface. The reference arm length is identical zero, so the cavity is build up by the gap between the Fizeau and the test surface only and contains no further optical elements. That is the reason why a Fizeau interferogram commonly directly shows the deviations of the test sample from the reference surface, i.e. Fizeau surface and allows skilled people analysis of the fringe pattern to judge the sample. The quality of the Fizeau surface determines the accuracy of the Fizeau interferometer. Fizeau surfaces are commonly available with a quality of  $\lambda/10$  -  $\lambda/20$  PV on special request also better.



## TRIOPTICS worldwide

### Headquarters

- TRIOPTICS Wedel - Germany

### Subsidiaries

- TRIOPTICS Berlin
- TRIOPTICS China
- TRIOPTICS France
- TRIOPTICS Hongkong
- TRIOPTICS Japan
- TRIOPTICS Korea
- TRIOPTICS Singapore
- TRIOPTICS Taiwan
- TRIOPTICS USA
- TRIOPTICS Wetzlar

### Partners

- Armstrong Optical - UK
- HP Instruments - India
- JSC Uran - Russia
- Optomek - Turkey
- Prolog Optics Ltd - Israel
- TECOTEC - Vietnam



TRIOPTICS Berlin GmbH • Interferometer Business  
Schwarzschildstraße 12 • 12489 Berlin / Germany  
Phone: +49-30-6392 3456  
Fax: +49-30-6392 3452  
E-mail: [support@trioptics-berlin.com](mailto:support@trioptics-berlin.com)  
[www.trioptics.berlin](http://www.trioptics.berlin)

2021 TRIOPTICS Berlin GmbH • All rights reserved

A member of the JENOPTIK Group