

## μPhase® VERTICAL

The universal interferometer system for production, workshop and laboratory.



### General features:

- Vertical setup with small footprint
- For measuring flat and spherical samples up to  $\varnothing = 50\text{mm}$
- Integrated automatic radius measurement
- Sample support on tilt and X-Y translation table
- Optional with tray for measuring a multiple number of samples at once
- Transmission measurement possible

If you need more information don't hesitate to contact us at [support@trioptics-berlin.com](mailto:support@trioptics-berlin.com).

(Picture above shows the version with μLens PLANO 50 for a single sample measurement.)