

Technical data ImageMaster® ENDO PRO

	ImageMaster® ENDO PRO
Accuracy (MTF on-axis and off-axis)	±0.03 MTF
Repeatability (MTF on-axis and off-axis))	±0.02 MTF
Endoscope lengths	140 mm ... 550 mm
Max. object distance	50 mm
Measurement time for all standard parameters	3 min
Measurement time for standard and extended Parameters	4 min
Off-axis measurement points	12 (others on request)
Direction of view	0 ... 70°
Max. field of view	120°
Max. object height	±35 mm
Spatial frequency	up to 350 lp/mm
Dimensions with table (height x width x depth)	1,900 mm x 1,550 mm x 900 mm
Weight with table	470 kg
Type	Stand alone