


μ Phase[®] ST / μ Phase[®] ST+R

A flexible and cost-effective solution for production and laboratory



Key Features:

- For measuring flat and spherical surfaces
- Available both for 10mm and 50mm based optics
- Table-top setup with small footprint
- Vertical column with manual coarse and fine z-movement
- Sample holder with tilt and X-Y translation table
- Optionally with radius measurement capability

Technical Data:

- Accuracy of radius measurement (optionally):
 - 0.5 μ m or 0.1 μ m (Heidenhain gauge, 60mm range)
 - 5 μ m (Mitutoyo gauge, 50mm range)
- Length of vertical column: 490 mm
- Footprint 280 mm x 360 mm
- Weight: ca 18 kg

For further information please contact support@trioptics-berlin.com.