

μPhase® PLANO DOWN



General features:

- Compact setup with small footprint
- For measuring flat samples up to 150mm in diameter
- Stable granite base in 2 sizes
- Tilt alignment integrated in the stand columns (micro-meter precision)

Available versions:

PLANO lens:	with 50mm / 100mm / 150mm free aperture Various adapters to exchange the lens in same stand
Stand size:	Small: Base 300 x 300 mm / free working width 180 mm Large: Base 400 x 400 mm / free working width 330 mm
Working height:	Selectable stand columns for free working height of 110mm / 155mm / 200mm

If you need more information don't hesitate to contact us at support@trioptics-berlin.com.

(Picture above shows the small version with 100mm μLens PLANO and 110 mm working height.)