

## μPhase® PLANO UP / μPhase® SPHERO UP

The μPhase® PLANO UP / μPhase® SPHERO UP are extremely compact and are ideally suited for production.



### General features:

- For measuring flat or spherical samples
- Small footprint
- Vibration insensitivity enables position directly next to the production machine
- Intuitive and easy handling enables usage by untrained personnel
- Suitable for integration into automated production lines

### Technical Data:

- μPhase® PLANO UP: flat surfaces  $\varnothing \leq 2\text{mm}$  to  $100\text{mm}$
- μPhase® SPHERO UP: convex from  $2\text{ mm}$  to  $225\text{ mm}$  radius and up to  $\varnothing = 55\text{ mm}$   
concave from  $-3\text{ mm}$  to  $-570\text{ mm}$  radius

If you need more information don't hesitate to contact us at [support@trioptics-berlin.com](mailto:support@trioptics-berlin.com).