

# μPhase® PLANO DOWN



## **Key Features:**

- Compact measurement setup with small footprint
- Stable granite base in two sizes
- For measuring flat surfaces up to 150mm in diameter
- Tilt alignment integrated in the stand columns (micro meter precision)

#### Available versions:

#### **PLANO** lens

- With 50mm / 100mm / 150mm free aperture
- Various adapters to change the lens in the same stand

### Stand size

- Small: Base 300 x 300mm / free working width 180mm
- Large: Base 400 x 400mm / free working width 330mm

## Working height

Selectable stand columns for free working height of 110 / 155 / 200mm

For further information please contact support@trioptics-berlin.com.