

Technical data ImageMaster® Afocal

	ImageMaster® Afocal
Spectral Range	450 ... 750 nm
Clear Aperture	25 ... 78 mm
Illumination Sources	Halogen light source (Optional "White" LED panel)
Max. off-axis object angle	±15 °
Max. off-axis image angle	±45° (90° optional)
Magnification	1 ... 60×
Accuracy (MTF on-axis and off-axis)	±3 MTF
Repeatability (MTF on-axis and off-axis)	±2 MTF
Object distances	20 m ... ∞
Dimensions (h x w x d)	253 mm x 1,782 mm x 627 mm
Weight	approx. 185 kg
Type	Table top device