

Technical data ImageMaster® MF

	ImageMaster® MF 100	ImageMaster® MF 200
Accuracy (MTF on-axis and off-axis)	±2 % MTF on-axis ±3 % MTF off-axis	±2 % MTF on-axis ±3 % MTF off-axis
Repeatability (MTF on-axis and off-axis)	±1 % MTF on-axis ±1,5 % MTF off-axis	±1 % MTF on-axis ±1,5 % MTF off-axis
Optical set up	Infinite conjugate set up	Infinite imaging set up
Sample EFL (Effective Focal Length)	1 mm ... 12 mm	18 mm ... 90 mm
Sample HFOV	10° ... 80°	6° ... 45°
Sample FOV	20° ... 160°	12° ... 90°
Dimensions (h x w x d)	1,570 mm x 1,340 mm x 810 mm	
Weight	240 kg	
Type	Stand alone	Stand alone